

# the Chair for Wound Care

The Hill 90W.

A dramatic upgrade from awkward stretchers, gurneys and chairs designed for other purposes.



*Shown with  
Azure upholstery  
and Beige base.*

## 90W Standard Features

- Height electrically controlled from 22"- 32"; Lifting capacity up to 600 lbs.
- Electric back with face cut-out and insert
- Electric tilt for Trendelenburg
- Electric foot section with 90° range
- Extending foot section - 8" range
- Width - 26"; Length - 71" (extends to 79")
- Comfortable contour back and seat with 1" soft foam over 1 1/2" medium density
- Durable upholstery; choice of 20 colors
- Baked on powder-coated base in 3 colors
- Met Lab/UL certified

## 90W Features and Options

**Hospital-grade Luxo Lamp**  
with magnifier and moveable stand



**Electric Back, Electric Foot and Electric Tilt**

The 90W easily converts from a seated position to Trendelenburg to Full recline.



▼ **Contoured back and seat**

▼ **Fold-up arm rests**

▼ **Matching Exam Stool**  
with or without back



▶ **Optional 30" wide** with removable armrest inserts



▶ **Breathing Slot** for Face-Down Positioning



▶ **Optl. Stainless Steel Tray**  
**Optl. Adjustable Stirrup** for Wrapping



**Quality Treatment Tables**  
Since 1945

  
**Hill Laboratories**  
COMPANY

Visit  
[www.HillLabs.com](http://www.HillLabs.com)  
for pricing and a video demonstration

1-877-445-5020 • [HillLabs.com/90w](http://HillLabs.com/90w)

Standard upholstery colors below with choice of base in beige, gray or black. Vinyl samples available upon request. For upgraded material and more information visit [HillLabs.com](http://HillLabs.com).

