

For Aesthetics

HA90ALB

The 90ALB (Aesthetic Lift Back) is an extremely affordable chair for skin care with 2" of soft foam over 1 1/2" of medium-density for maximum patient comfort. Many standard and optional features make the 90ALB ideal for aesthetic work.



- Shown with
Bordeaux upholstery
and Beige base.

Easy up. Easy down.

The 90ALB features electric elevation from 23" to 36" or 21" to 31". The standard model has an assisted, gas spring back for easy positioning, an adjustable headrest, fold-down movable armrests, contour back and durable, upgraded designer upholstery. The HA90ALB can fully recline and is also available with an optional electric lift back and a variety of head cushion options.

Powered by an electric actuator 60 Hz,
110 volts. Amperage - 5.1

Lifting Capacity:

400 lbs. (600 lbs. capacity upon request)

Height Range:

90ALB - 23" to 36" or 21" to 31"

Table Base:

Beige, Gray or Black telescoping base.

All metal parts are powder coated.

Shipping Weight:

90ALB - 225 lbs.

Warranty: One year



Hill Laboratories
COMPANY

1-877-445-5020 • www.HillLabs.com