



Hill Massage Table



Comfortable for you and your clients with a safer, stronger design.


Just \$1895!

A PROFESSIONAL LOOK at an INCREDIBLE PRICE!

The attractive design of the Hill Massage Table inspires confidence and breathes professionalism. It's fully enclosed base is safer and simply looks better than many of the open-base tables on the market today.

Hill Massage Table Features

- Electric Power Height from 21" to 31"
- Height controlled by foot pedal
- Custom Comfort Top - 27" or 30" wide with 1" soft over 1.5" medium density foam combination.
- Easy-Lock Face Piece - locks at any angle positive or negative
- Plush and durable vinyl upholstery available in 8 colors with PreFixx® stain resistance
- Fully enclosed base available in beige, gray or black
- One-year warranty on all parts (Dealer on-site labor warranty may vary)



Available Colors for the Hill Massage Table.



Hill Laboratories
COMPANY

Call now! 1-877-445-5020
or visit www.HillLabs.com