

# Hill Laboratories

## MAT TABLE

### Electric Elevation HA90M



### HA90M FEATURES

Hill Laboratories has designed the HA90M to be sturdy, dependable, affordable—and the most attractive mat table on the market! All of its working mechanisms are safely enclosed within an acrylic, thermo-plastic base which is durable and very easy-to-clean. The base features a welded steel construction and a baked-on, powder-coat finish. Vinyl bumper guards protect the edges of the table-top while recessed support legs permit obstruction-free movement around the table. The Hill Adjustable Mat Table adjusts smoothly and quietly from 21" to 34" with the power of an electric actuator. Design the HA90M with your choice of 20 attractive Blockade® upholstery colors.

- 2" medium-density Ultra Cell foam for comfort & support
- Standard top sizes of 4' x 6', 4' x 7', 5' x 7', 6' x 8
- Custom top sizes also available
- Standard height of 21" - 34"
- Special low height of 19" available
- Lifting capacity of 700 lbs.
- Elevates from low to high in 15 seconds
- Choice of gray or beige base
- Vinyl bumper guard surrounds perimeter of table top
- Upholstered with highly-durable, low-maintenance Blockade material; also available in exceptionally strong Herculite
- Optional extra foot pedal
- Electric requirement- 110 volt, 60 Hz, 2.7 amp.
- One-year warranty on parts and labor
- UL Listed

  
Hill Laboratories  
COMPANY