

**Saves Time. Better Treatments.**  
Multi-position soundhead for consistent, repeatable hands-free treatments.

**Versatile**

Includes interferential and premod current; use alone or in combination with ultrasound directly through the soundhead.



**Powerful Infrared.**

Optional as a plug-in



**Optional Handheld Applicator**



**Accessories.**

Comes complete with 2 Muscle Stimulation Wires, Electrodes, Gel Pads, Weight Bag, Velcro Straps.



**Mobile.**

Rolling Cart with Basket included.

# HF54 PLUS

- ▶ Hands-Free & Optional Hand-held Ultrasound
- ▶ Premod, Interferential, Russian, & TENS
- ▶ Combo Mode directly through the Soundhead
- ▶ Add Anytime: Optional 2nd Soundhead and Light Therapy

The HF54 Plus features three uniform harmonized crystals that produce an even distribution of sound waves throughout its 3-5/8 diameter. The HF54 operates at 1 MHz; the crystals pulsate at 50% ON/50% OFF at 60 Hz.

Following our guide for intensity enables the user to leave the sound head safely on an area for the entirety of the treatment time.

### Combination Therapy

Interferential Russian or TENS can be used in one area while treating with Ultrasound in another.

Ultrasound can be used in combination with premod allowing for ultrasound and muscle stim through soundhead.

The reusable gel pad is the conductive material for the soundhead making for easy patient setup.

- ▶ *See other side for details*



Hill Laboratories  
COMPANY

1-877-445-5020  
www.HillLabs.com

# HF54

# A Unique Soundhead

## Standard Features

- HF54 PLUS Hands-Free Therapeutic Ultrasound Generator
- Swing-Arm with Soundhead
- IFC, Premod Current, TENS and Russian
- 2 Muscle Stimulation Wires
- Two Packs of 2.75" x 4.75" Electrodes and Two Packs of 2" x 2" Electrodes
- Cut Off Switch
- 4 Gel Pads & 1 Retainer Ring
- Weight Bag & Velcro Straps
- Rolling Cart with Basket
- Utility Tray & Spray Bottle
- FDA Cleared
- One-year warranty on the control unit and soundhead

Easily add a second soundhead at any time for treating two areas at the same time



▲ Back - Ultrasound with Dual Soundheads

## Apply with Weight Bag or specially-designed Velcro Straps



▲ Ultrasound on Cervical Region



▲ Ultrasound for Foot



▲ Ultrasound for Elbow



▲ Ultrasound with Premod Combination



▲ Ultrasound for Knee



▲ Interferential on Lower Back



## Hand-Held Ultrasound

Add optional hand-held ultrasound to the HF54 Plus at any time with a traditional 5 cm. soundhead which can be used at 1 or 3 MHz.; continuous or pulsating.



## Horizon Light Therapy Unit

Add optional light therapy to your HF54 PLUS. The Horizon features a two inch diameter treatment area, 800mW Infrared (880nm) and Red (660nm) light with 36 infrared super luminous diodes and 7 visible red super luminous diodes.



See the video at [HillLabs.com](http://HillLabs.com)

  
**Hill Laboratories**  
COMPANY

1-877-445-5020  
[www.HillLabs.com](http://www.HillLabs.com)