


Hill Adjustable

# DNFT TABLE

***Directional Non-Force Technique***



*Pictured in Azure upholstery  
with a gray base.*

## Featuring

- Variable height from 22" to 35"
- Auto Touch foot pedal
- Tilting headpiece
- 13" extended Face-Cut out into Thoracic Cushion with Insert Filler
- Slide-out ankle support extends 12"; adjusts in height 8"
- Measurements - 5'7" - 6'7" long; 24" wide
- Arm rests
- Paper-roll bracket
- 2 1/2" medium-density foam top; head cushion with soft foam
- Choice of 20 upholstery colors
- Beige, gray or black acrylic base

Equipped with an Autotouch foot pedal which will raise the table to the maximum height AND lower the table to the minimum height respectively when tapped a single time. Operators may stop the table at any point within the travel cycle.

"D.N.F.T.® is the exclusive registered trademark of Dr. Chris John, teacher of Directional Non-Force Technique® Chiropractic. The "Hill D.N.F.T. Elevation Table" is developed through the collaboration of Hill Laboratories and Dr. Chris John, and manufactured by Hill Laboratories with Dr. John's design approval."

Freight charges are additional • Prices are subject to change without notice.

  
**Hill Laboratories**  
**COMPANY**

3 Bacton Hill Road, Frazer, PA 19355 • 1-877-445-5020 • 610-644-2867 • fax 610-647-6297 • [www.HillLabs.com](http://www.HillLabs.com)