

# Electric Adjustable Height

## HILL CAST TABLE

### PERFECT PATIENT POSITIONING AND TRANSFER



The Hill Adjustable Cast Table (90CT) is designed to make your work easier.

With elevation from 22" to 35" and rocker foot pedals to adjust the height from either side of the table, the Hill Cast Table offers easy patient transfer and the ideal working height for you. A gas-spring assisted lift-back raises with ease to any position up to 80°

The leg section on the 90CT is split offering independent right and left leg adjustment at any angle down to 90°. A foot plate greatly simplifies short leg casts and posterior ankle splints.

The Hill Cast Table features fold-down armrests which also swing up and out of the way when not in use. The armrests have a removable rubber handle for thumb spica splints and other finger, thumb and wrist splints. An optional rotating armboard which moves and locks on both the X and Y axis, can also be added. The Hill Cast Table has a one-year warranty and is available with a gray or beige base and several attractive upholstery colors.



**Hill Laboratories**  
COMPANY

1-877-445-5020 • [www.HillLabs.com](http://www.HillLabs.com)