

the Chair for Wound Care

The Hill 90W.

A dramatic upgrade from awkward stretchers, gurneys and chairs designed for other purposes.



Shown with
Azure upholstery
and Beige base.

90W Standard Features

- Height electrically controlled from 22" - 32"; Lifting capacity up to 600 lbs.
- Electric back with face cut-out and insert
- Electric tilt for Trendelenburg
- Electric foot section with 90° range
- Extending foot section - 8" range
- Width - 26"; Length - 71" (extends to 79")
- Comfortable contour back and seat with 1" soft foam over 1 1/2" medium density
- Durable upholstery; choice of 20 colors
- Baked on powder-coated base in 3 colors
- Met Lab/UL certified

90W Features and Options

Hospital-grade Luxo Lamp
with magnifier and moveable stand



Electric Back, Electric Foot and Electric Tilt
The 90W easily converts from a seated position to Trendelenburg to Full recline.



- ▼ Contoured back and seat
- ▼ Fold-up arm rests
- ▼ Matching Exam Stool with or without back



Optional 30" wide with removable armrest inserts



Breathing Slot for Face-Down Positioning



Optl. Stainless Steel Tray
Optl. Adjustable Stirrup for Wrapping



Quality Treatment Tables
Since 1945


Hill Laboratories
COMPANY

Visit
www.HillLabs.com
for pricing and a video demonstration

1-877-445-5020 • HillLabs.com/90w

Standard upholstery colors below with choice of base in beige, gray or black. Vinyl samples available upon request. For upgraded material and more information visit HillLabs.com.

