

The Counterstrain Table

by Hill Laboratories



Featuring-

- Tilting Headpiece - 90° negative and 30° positive
- Electric Tilt-Back
- Slide-out pelvic section
- Lift leg section - gas spring assisted
- Drop-foot section, 8", Mechlok® operated, 90° negative and 30° positive.
- Double locks (Mechlok) on each side of the table to raise and lower head and leg section
- Two foot pedals mounted near the front of the table to control head/thoracic area
- Two foot pedal mounted near the base of the table to control height
- Headpiece open at forehead area
- Armrests mounted on head section
- Electric hand-operated switch to elevate table at head end
- 78" long, 24" wide
- 20 designer colors available for the top and gray, beige and black for the base


Hill Laboratories
COMPANY